




# Product Briefing

## Chromacol 9-SCJ caps

Caps and seals for Shimadzu and Jasco Instruments

- Fits the wide range of Chromacol wide opening screw thread vials.
- The cap has a support ring to allow its use in standard 12mm vial racks in early and late model Shimadzu SIL-6 and SIL-10 auto samplers and also Jasco instruments.



Chromacol 9-SCJ caps combine the auto sampler capability of the industry standard 9mm vials with Japanese instruments.

Chromacol 9-SCJ caps are designed to fit the full range of Chromacol 9mm vials from the standard 2mL to High Recovery and Micro+™ vials.

This allows use of the popular 03-FISV vials.

### ***Chromacol 9-SCJ Caps provide:-***

- Broad compatibility with robotic auto sampler vials including Micro+™.
- Pre-inserted Rubber/PTFE and Silicone/PTFE seals.
- Specific colouring to give the correct vial detection in photo diode identification systems.

### ***Chromacol 9-SCJ Caps***

Chromacol Cap	Colour	Cap Material	Vial Compatibility
9-SCJ	Black	HDPE	9mm Screw
9-SCJ(Y)	Yellow	HDPE	9mm Screw
9-SCJ(W)	White	HDPE	9mm Screw

Chromacol 9mm Screw Vial	Description
2-SVW	2mL wide mouth screw top vial with write on patch - clear
2-SVW(A)	2mL wide mouth screw top vial with write on patch - amber
1.5-HRSV	1.5mL screw top vial - clear - High Recovery **Available in packs of 100**
1.5-HRSV(A)	1.5mL screw top vial - amber - High Recovery **Available in packs of 100**
09-FISV	.9mL screw top fixed insert vial - clear ( Micro+™)
03-FISV	.3mL screw top fixed insert vial - clear with write on patch ( Micro+™)
03-FISV(A)	.3mL screw top vial fixed insert - amber with write on patch ( Micro+™)



# Product Briefing

## Chromacol 9-SCJ caps and seals

Chromacol Cap	Colour	Description	Compatibility
9-SCJ-8RT1	Black	9mm screw cap with flange for Japanese instruments pre-fitted type 8 rubber/PTFE liner	General purpose use with 2-SVW and other 9mm robotic type vials
9-SCJ-ST15	Black	9mm screw cap with flange for Japanese instruments pre-fitted silicone/PTFE liner	General purpose use with 2-SVW and other 9mm robotic type vials.
9-SCJ-ST101	Black	9mm screw cap with flange for Japanese instruments pre-fitted soft silicone/PTFE liner	High purity soft seal for use with MS and Specific detectors
9-SCJ(W)-8RT1	White	9mm white screw cap with flange for Japanese instruments pre-fitted type 8 rubber/PTFE liner	Shimadzu LC Autosamplers SIL-10, SIL-6 General purpose use with 2-SVW and other 9mm robotic type vials
9-SCJ(W)-ST15	White	9mm white screw cap with flange for Japanese instruments pre-fitted silicone/PTFE liner	Shimadzu LC Auto samplers SIL-10, SIL-6 General Purpose use with 2-SVW and other 9mm robotic type vials
9-SCJ(W)-ST101	White	9mm white screw cap with flange for Japanese instruments pre-fitted soft silicone/PTFE liner	Shimadzu LC Auto samplers SIL-10, SIL-6 High purity seal for use with 2-SVW and other 9mm robotic type vials
9-SCJ(Y)-8RT1	Yellow	9mm yellow screw cap with flange for Japanese instruments pre-fitted type 8 rubber/PTFE liner	Jasco AS-1500/2000 series General purpose use with 2-SVW and other 9mm robotic type vials
9-SCJ(Y)-ST15	Yellow	9mm yellow screw cap with flange for Japanese instruments pre-fitted silicone/PTFE liner	Jasco AS-1500/2000 series General purpose use with 2-SVW and other 9mm robotic type vials
9-SCJ(Y)-ST101	Yellow	9mm yellow screw cap with flange for Japanese instruments pre-fitted soft silicone/PTFE liner	Jasco AS-1500/2000 series High purity seal for use with 2-SVW and other 9mm robotic type vials

Chromacol Ltd  
3 Mundells Industrial Centre  
Welwyn Garden City  
Herts AL7 1EW  
United Kingdom

Phone +44 (0) 1707 394949  
enquiries@chromacol.com

Fax +44 (0) 1707 391311  
www.chromacol.com