



# Chromacol Headspace Vial Selection



Manufacturer	Instrument Ranges	Cap Parts	Vials		Instrument	
Agilent	7694 Headspace	 20-AC-ST3(W)	 20-CV	 9-CV		
Agilent	G1888	 20-AC-ST3(W)	 20-CV	 10-CV		
CTC	Combi PAL	 20-MCBC-ST3	 20-CV	 10-CV		
DANI	HSS 39.50	 20-AC-CBT3	 20-CV	 9-CV		
DANI	HSS 86.50	 20-AC-CBT3	 20-CV	 9-CV		
NLG/ Eurovector/ SRI	HT200H	 20-AC-ST3	 20-CV	 10-CV		
PE	HS16	 20-AC-ST3	 20-CV	 10-CV		
PE	HS40	 20-AC-ST3	 20-CV	 6-CV		
PE	HS110	 20-AC-ST3	 20-CV	 6-CV		

# Chromacol Headspace Vial Selection



Manufacturer	Instrument Ranges	Cap Parts	Vials			Instrument	
Shimadzu	HSS-2B	 20-AC-ST3	 27-CV				
Shimadzu	AOC-5000	 20-MCBC-ST3	 20-CV	 10-CV			
Teledyne Tekmar	7000/7050	 20-AC-CBT3	 22-CV	 12-CV	 9-CV		
Thermo/ Fisons/ Carlo Erba	HS500/ HS850	 20-MCB and 20-ST15	 20-CV (HS850)	 9-CV (HS500)			
Thermo	HS2000	 20-AC-ST3	 20-CV	 10-CV			
Thermo	Finnigan Triplus	 20-MCBC-ST3	 20-CV	 10-CV			
Varian	Genesis	 20-AC-CBT3	 22-CV	 12-CV	 9-CV		
Varian	CP8000 Series Combi PAL	 20-MCBC-ST3	 20-CV	 10-CV			