



Product Briefing

Chromacol Magnetic Cap Adapter 9-MC

New Accessory

With new sample processors and autosamplers being developed to move vials using magnetic attraction one of the requests is a means of transporting 2mL vials to allow their use in SPME and headspace extraction methods. Until now this has required the use of tinplate caps. Now a simple reusable accessory can convert 2mL vials into a transport vial. The magnetic component is a push-fit 9-MC that is re-usable, allowing repeated use of the transfer facility. The vial is a 2-SVW and can be combined with the 9mm flanged screw cap, 9-SCJ. These can be combined with the high purity ST-101 seal.

The external dimensions allow the use of these in trays and racks designed for a number of units including the Thermo Finnigan Triplus HS, the CTC CombiPAL and Gerstel MPS-2.



Ordering

Part	Description
9-MC	Magnetic sleeve for 9-SCJ Cap-Pk-100
9-SCJM-ST101	9mm Magnetic Screw Cap Assembly with Pre-fitted Soft Silicone/PTFE Seal Pk-100
9-SCJ-ST101	9mm Screw Cap with Flange with Pre-fitted Soft Silicone/PTFE Seal
2-SVW	2mL Wide mouth screw top vial with write on patch – clear

Chromacol Ltd
3 Mundells Industrial Centre
Welwyn Garden City
Herts. AL7 1EW
United Kingdom

Phone +44 (0) 1707 394949
Fax +44 (0) 1707 391311

Email
enquiries@chromacol.com

Website
www.chromacol.com