

Chromacol Instrument SELECT Packs

GC



CHROMACOL




Part	Instrument Ranges	Components	Chromacol Parts	Vial	Instrument	
CTCG	CTC/Leap GCPAL CombiPAL	- 2mL crimp top vial with a blue silicone/red PTFE seal	2-CV 11-AC-ST144	 2-CV		
HPG	Agilent Hewlett Packard	- 2mL crimp top vial with a type 7 rubber/PTFE seal	2-CV 11-AC7	 2-CV		
PEG	Perkin Elmer AutoSystem	2mL crimp top vial with a type 6 soft rubber/PTFE seal.	2-CV 11-AC6	 2-CV		
SHG	Shimadzu	- 2mL crimp top vial with a soft blue silicone/PTFE seal	2-CV 11-AC-ST101	 2-CV		
TQG	Thermo Electron Thermo Finnigan ThermoQuest TSP Carlo Erba Fisons Instruments	- 2mL crimp top vial with a blue silicone/red PTFE seal	2-CV 11-AC-ST144	 2-CV		
TTR	Thermo Electron Thermo Finnigan TRACE	2mL screw top vial with a blue screw cap with soft blue silicone/PTFE seal	2-SVW 9-SC(B)-ST101	 2-SVW		
VAG	Varian Chrompack Packard	- 2mL screw top vial with a white silicone/red PTFE seal	2-SVW 9-SC(B)-ST1	 2-SVW		

Chromacol SELECT Packs

IC



CHROMACOL

Part	Instrument Ranges	Components	Chromacol Parts	Vial	Instrument	
DIC	Dionex Ion Chromatography AS40/AS50	0.6mL screw top polypropylene vial with pre-slit silicone/PTFE seal	06-PESV 8-SC-ST14X	 06-PESV		

Chromacol Instrument SELECT Packs

HPLC



Part	Instrument Ranges	Components	Chromacol Parts	Vial	Instrument	
CTCL	CTC/Leap LC PAL	- 2mL crimp top vial with a blue cap with white silicone/red PTFE seal	2-CV 11-AC(B)-ST144	 2-CV		
HPL	Agilent 1100 Hewlett Packard 1050 1090	- 2mL crimp top vial with a type 7 rubber/PTFE seal	2-CV 11-AC7	 2-CV		
HPLS	Agilent 1100 Hewlett Packard 1050 1090	2mL screw top vial with a type 8 rubber/PTFE seal	2-SVW 9-SC(B)-8RT1	 2-SVW		
MEL	Merck-Hitachi LaChrom	- 2mL crimp top vial with a pre-cut blue silicone/PTFE seal	2-CV 11-AC-ST101X	 2-CV		
PEL	Perkin Elmer Series 200	- 2mL screw top vial with a pre-cut white silicone/red PTFE seal	2-SVW 9-SC(B)-ST1X	 2-SVW		
SHL	Shimadzu	- 2mL screw top vial with a soft blue silicone/PTFE seal and a wider screw cap	2-SV 8-SCJ(W) 8-ST101	 2-SV		
SPL	Spark	- 2mL crimp top vial with a type 7 rubber/PTFE seal	2-CV 11-AC7	 2-CV		
TQL	Thermo Finnigan ThermoQuest TSP Spectraphysics LDC	- 2mL crimp top vial with a white silicone/red PTFE seal	2-CV 11-AC-ST15	 2-CV		
VAL	Varian	- 2mL screw top vial with a white silicone/red PTFE seal	2-SVW 9-SC(B)-ST1	 2-SVW		
WAL	Waters Alliance	- 2mL screw top vial with a white silicone/red PTFE seal	2-SVW 9-SC(B)-ST1	 2-SVW		

All Trademark rights are acknowledged