

Agilent Technologies (HP) 1050, 1100 and 1200

Our 11-AC7, red rubber and PTFE, crimp caps, are ideal for this instrument and our 11-AC-ST15 and 11-AC-ST101 silicone seals may be used when a more aggressive solvent is specified.

The Instrument Select product code HPL includes 1000 2-CV vials and 1000 11-AC7 crimp caps specially selected to be compatible with this instrument.

Key to products ■ Vial ■ Cap ■ Combination Pack

Bold items in grey are our standard products, used in the combination packs.

| PART NUMBER | DESCRIPTION | 12 x 32mm crimp and snap cap vials and caps | PACK SIZE |
|-----------------------|--|---|------------|
| HPL | 2mL crimp top vial with a type 7 rubber/PTFE seal. | | 100 |
| 2-CV | 2mL, clear glass, crimp top, flat bottom, wide neck vial. With write-on patch. | | 500 |
| 1.1-CTVG | 1.1mL, clear glass, crimp top, tapered, wide neck vial. Made from inert Chromacol Gold™ grade glass. Needs PTFE support, reference TTS-313. See page 78. | | |
| 09-FIV | 900µL Micro+ crimp top, fused insert vial, clear. | | 500 |
| 06-PECV | 600µL, polyethylene vial with an internal tapered profile. | | 500 |
| 03-FIV | 300µL, Micro+ one piece, clear glass insert and vial combination. With write-on patch. | | 500 |
| 11-AC7 | Aluminium crimp cap fitted with a type 7, rubber/PTFE seal. Shore hardness value 60. Other colours and seals available. See page 32. | | 500 |
| 2-RV | 2mL snap cap vial with write on patch - clear. | | 500 |
| 11-PSN(B)-8RTI | 11mm polyethylene snap cap pre-fitted with a red rubber/PTFE seal. Shore hardness value 58. | | 500 |

| PART NUMBER | DESCRIPTION | 6 and 7 x 32mm crimp top vials and cap | PACK SIZE |
|------------------|---|--|-----------|
| 08-CRV(A) | 800µL, amber glass, crimp top, round bottom vial. Needs WS-5 support. | | 500 |
| 06-CTV(A) | 600µL, amber glass, crimp top, tapered vial. Needs WS-5 support. | | 500 |
| 03-CVG | 300µL, clear glass, crimp top, round bottom vial, made from inert Chromacol Gold grade glass. Needs support reference SV-S11A. See page 28. | | |
| 8-AC7 | Aluminium crimp cap fitted with a type 7, rubber/PTFE seal. Shore hardness value 60. Other colours and seals available. See page 29. | | 1000 |

| PART NUMBER | DESCRIPTION | 12 x 32mm wide neck screw top vials and cap | PACK SIZE |
|---------------------|--|---|-----------|
| HPLS | 2mL screw top vial with a type 8 rubber/PTFE seal. | | 100 |
| 2-SVW | 2mL, clear glass, flat bottom, screw top, wide neck vial. With write-on patch. | | 500 |
| 09-FISV | 900µL Micro+, screw top, fused insert vial, clear. | | 500 |
| 03-FISV | 300µL Micro+, one piece, clear glass insert and vial combination with write-on patch. | | 500 |
| 9-SC(B)-8RTI | 9mm screw cap (blue) prefitted rubber PTFE seal. Shore hardness value 58. Other colours and seals available. See page 36. | | 500 |
| 4-SV | 4mL wash vial, 15 x 46mm, clear glass, screw top, flat bottom. See page 38. | | 500 |

