

Agilent Technologies (HP) 7673A/1/1/1

One of the world's most popular autosamplers is the HP 7673A, and there are many Chromacol vials of different shape and volume, caps and seals with many different formulations, which are ideal for this instrument. Our 2-CV vials and 11-AC7 caps, included in our Instrument Select packs HPG, are a perfect match and our Micro+ range of vials are ideal for microsampling. For the full range of our GC Septa, please request a copy of our GC Septa brochure or see page 56.

Key to products ■ Vial ■ Cap ■ Combination Pack

Bold items in grey are our standard products, used in the combination packs.

PART NUMBER	DESCRIPTION	12 x 32mm crimp and snap cap vials and caps	PACK SIZE
HPG	2mL crimp top vial with a type 7 rubber/PTFE seal.		100
2-CV	2mL, clear glass, crimp top, flat bottom, wide neck vial. With write-on patch.		500
09-FIV	900µL Micro+ crimp top, fused insert vial, clear.		500
06-PECV	600µL, polyethylene vial with an internal tapered profile.		500
03-FIV	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
11-AC7	Aluminium crimp cap fitted with a type 7, rubber/PTFE seal. Shore hardness value 60. Other colours and seals available. See page 32.		500
2-RV	2mL snap cap vial with write on patch - clear.		500
11-PSN(B)-8RT1	11mm polyethylene snap cap pre-fitted with a red rubber/PTFE seal. Shore hardness value 58.		500

PART NUMBER	DESCRIPTION	6 x 32mm crimp top vial and cap	PACK SIZE
03-CVG	300µL, clear glass, crimp top, round bottom vial, made from inert Chromacol Gold grade glass, needs SV-S11G or SV-S11A support. See page 28.		500
8-AC7	Aluminium crimp cap fitted with a type 7, rubber/PTFE seal. Shore hardness value 60. Other colours and seals available. See page 29.		1000

PART NUMBER	DESCRIPTION	12 x 32mm wide neck screw top vials and cap	PACK SIZE
2-SVW	2mL, clear glass, flat bottom, wide neck, screw top vial. With write-on patch.		500
09-FIV	900µL Micro+ crimp top, fused insert vial, clear.		500
03-FISV	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
9-SC(B)-8RT1	9mm screw cap (blue) prefitted rubber PTFE seal. Shore hardness value 58. Other colours and seals available. See page 36.		500
4-SV	4mL wash vial, 15 x 46mm, clear glass, screw top, flat bottom. See page 38.		500

Chromaseal™ GC Septa	Model	Septa size		PART NUMBER			Pack Size
		mm	ins.	High Temp	Long Life	Economy	
Agilent Technologies (HP)	5700	9.5	3/8	HT-9.5	LL-9.5	ECO-9.5	25
	5880, 5890, 6890	11	7/16	HT-11	LL-11	ECO-11	25

