

Agilent Technologies (HP) 1090

Millions of Chromacol vials have been used with this instrument for many years in its characteristic black cassette system. Our Micro+ vials with a capacities of 300µL and 900µL are both easy to use and very safe for the autosampler needle.

Key to products ■ Vial ■ Cap ■ Combination Pack

Bold items in grey are our standard products, used in the combination packs.

PART NUMBER	DESCRIPTION	12 x 32mm crimp and snap cap vials and caps	PACK SIZE
HPL	2mL crimp top vial with a type 7 rubber/PTFE seal.		100
2-CV	2mL, clear glass, crimp top, flat bottom, wide neck vial. With write-on patch.		500
1.1-CTVG	1.1mL, clear glass, crimp top, tapered, wide neck vial. Made from inert Chromacol Gold grade glass. Needs PTFE support, reference TTS-313. See page 78.		500
09-FIV	900µL Micro+ crimp top, fused insert vial, clear.		500
06-PECV	600µL, polyethylene vial with an internal tapered profile.		500
03-FIV	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
11-AC7	Aluminium crimp cap fitted with a type 7, rubber/PTFE seal. Shore hardness value 60. Other colours and seals available. See page 32.		500
2-RV	2mL snap cap vial with write on patch - clear.		500
11-PSN(B)-8RT1	11mm polyethylene snap cap pre-fitted with a red rubber/PTFE seal. Shore hardness value 58.		500

PART NUMBER	DESCRIPTION	6 x 32mm crimp top vial and cap	PACK SIZE
03-CVG	300µL, clear glass, crimp top, round bottom vial, made from inert Chromacol Gold grade glass needs SV-S4 support. See page 28.		500
8-AC7	Aluminium crimp cap fitted with a type 7, rubber/PTFE seal. Shore hardness value 60. Other colours and seals available. See page 29.		1000

PART NUMBER	DESCRIPTION	12 x 32mm wide neck screw top vials and cap	PACK SIZE
2-SVW	2mL, clear glass, flat bottom, screw top, wide neck vial. With write-on patch.		500
09-FISV	900µL Micro+ screw top, fused insert vial, clear.		500
03-FISV	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
9-SC(B)-8RT1	9mm screw cap(blue) prefitted rubber PTFE seal. Shore hardness value 58. Other colours and seals available. See page 36.		500
4-SV	4mL wash vial, 15 x 46mm, clear glass, screw top, flat bottom. See page 38.		500

