

The 11-AC-ST101X and 8-ST14X products are suited to these autosamplers that require the use of a pre-cut seal. Merck/Hitachi have a wide variety of racks so they are compatible with a wide variety of vials. The 03-CVG, 02-CTVG, 02-CTV(A) and 01-CVG vials may be used in some trays without a support sleeve. The 2-CV crimp top vial and the 11-AC-ST101X caps are included in the Instrument Select kit – product code MEL.

**Key to products** ■ Vial ■ Cap ■ Seal ■ Combination Pack

**Bold items in grey are our standard products**

PART NUMBER	DESCRIPTION	15 x 46mm screw top vials and caps	PACK SIZE
<b>4-SV</b>	4mL, clear glass, screw top, flat bottom vial.		500
<b>12-SC-ST2</b>	12mm, black, screw cap pre-fitted with a silicone/PTFE seal. Shore hardness value 57.		500
<b>3.5-HRSV</b>	3.5mL screw top vial - High Recovery - uses 13mm screw caps.		500
<b>13-SC-ST15X</b>	13mm screw cap prefitted silicone/PTFE seal, pre-cut.		500

PART NUMBER	DESCRIPTION	12 x 32mm crimp and snap cap vials and caps	PACK SIZE
<b>MEL</b>	2mL crimp top vial with a pre-cut blue silicone/PTFE seal.		100
<b>4-CV</b>	4mL, clear glass, crimp top, flat bottom vial, 12 x 46mm.		500
<b>2-CV</b>	<b>2mL, clear glass, wide neck, crimp top, flat bottom vial. With write-on patch.</b>		<b>500</b>
<b>1.1-CTVG</b>	1.1mL, clear glass, crimp top, tapered vial made from inert Chromacol Gold grade glass, needs a PTFE support sleeve, part number TTS-312. <a href="#">See page 78.</a>		500
<b>09-FIV</b>	900µL Micro+, crimp top, fused insert vial, clear.		500
<b>06-PECV</b>	600µL Micro+, polyethylene, crimp top vial with an internal tapered profile.		500
<b>03-FIV</b>	300µL, one piece, clear glass insert and vial combination. With write-on patch.		500
<b>11-AC-ST101X</b>	<b>Aluminium crimp cap fitted with a pre-cut blue silicone/PTFE seal.</b>		<b>500</b>
<b>2-RV</b>	2mL snap cap vial with write on patch - clear.		500
<b>11-PSN(B)-STIX</b>	11mm polyethylene snap cap with pre-cut silicone/PTFE seal.		500

PART NUMBER	DESCRIPTION	12 x 32mm narrow & wide neck screw top vials and caps	PACK SIZE
<b>2-SVW</b>	2mL, clear glass, flat bottom, screw top, wide neck vial. With write-on patch.		500
<b>09-FISV</b>	900µL Micro+ screw top, fused insert vial, clear.		500
<b>03-FISV</b>	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
<b>9-SC(B)-STIX</b>	Screw cap, blue, pre-fitted with a blue silicone/PTFE pre-cut seal. Other colours and seals available. <a href="#">See page 36.</a>		500
<b>2-SV</b>	2mL screw top vial with write on patch - clear.		500
<b>1.1-STVG</b>	1.1mL, clear glass, 8mm, screw top, tapered vial made from inert Chromacol Gold grade glass, needs a TTS-312 support. <a href="#">See page 78.</a>		500
<b>8-SC</b>	8mm screw cap, black. For use with 2-SV or 1.1-STVG		500
<b>8-ST14X</b>	8mm pre-cut silicone/PTFE seal.		500

PART NUMBER	DESCRIPTION	6 x 32mm crimp top vial and cap	PACK SIZE
<b>03-CVG</b>	300µL, clear glass, crimp top, round bottom vial, made from inert Chromacol Gold grade glass, needs support sleeve SV-S1. <a href="#">See page 28.</a>		500
<b>8-AC-ST101X</b>	Aluminium crimp cap fitted with a pre-cut blue silicone/PTFE seal.		500

