

PerkinElmer AutoSystem/XL and Clarus 500/600

These autosamplers need to be used with relatively soft seals and the Chromacol type 6 rubber/PTFE seal was developed specifically for these instruments. Our silicone/PTFE formulation ST101 has been proven to be ideal for use with ECD on this instrument. Instrument Select pack PEG contains the standard caps and seals for this instrument. For the full range of our GC Septa, please request a copy of our GC Septa brochure or see page 56.

Key to products ■ Vial ■ Cap ■ Combination Pack

Bold items in grey are our standard products, used in the combination packs.

PART NUMBER	DESCRIPTION	12 x 32mm crimp top vials and caps	PACK SIZE
PEG	2mL crimp top vial with a type 6 soft rubber/PTFE seal.		100
2-CV	2mL, clear glass, crimp top, flat bottom, wide neck vial. With write-on patch.		500
09-FIV	900µL Micro+ crimp top, fused insert vial, clear.		500
06-PECV	600µL, polyethylene, crimp top vial with an internal tapered profile.		500
03-FIV	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
11-AC6	Aluminium crimp cap fitted with a type 6, rubber/PTFE seal. Shore hardness value 38. Other colours and seals available. See page 32.		500
11-AC-ST101	Aluminium crimp cap fitted with a blue silicone/PTFE seal. Shore hardness value 30, also good for ECD. Other colours and seals available. See page 32.		500

PART NUMBER	DESCRIPTION	6 x 32mm crimp top vial and cap	PACK SIZE
03-CVG	300µL, clear glass, crimp top, round bottom vial, made from inert Chromacol Gold grade glass. Needs a SV-TSP support. See page 28.		500
8-AC6	Aluminium crimp cap fitted with a type 6, rubber/ptfe seal. Shore hardness value 38. Other colours and seals available. See page 29.		1000

PART NUMBER	DESCRIPTION	12 x 32mm wide neck screw top vials and cap	PACK SIZE
2-SVW	2mL, clear glass, flat bottom, screw top, wide neck vial. With write-on patch.		500
09-FISV	900µL Micro+ screw top, fused insert vial, clear.		500
03-FISV	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
9-SC(B)-ST101	9mm, screw cap pre-fitted with a blue silicone/PTFE seal, Shore hardness value 30.		500

Chromacol GC Septa	Model	Septa size		PART NUMBER			Pack Size
		mm	ins.	High Temp	Long Life	Economy	
PerkinElmer	All models	11	7/16	HT-11	LL-11	ECO-11	25

