

PerkinElmer ISS-100, 200, Integral 4000 and Series 200

PerkinElmer LC autosamplers use a series of racks to hold the vials and there are many Chromacol vials compatible with them from 100µL to 10mL. Chromacol's WebSeal 96 well titer plates products are also suitable for certain versions of these instruments. To order easily, use Instrument Select code PEL.

Key to products ■ Vial ■ Cap ■ Combination Pack

Bold items in grey are our standard products, used in the combination packs.

PART NUMBER	DESCRIPTION	12 x 32mm crimp and snap cap vials and caps	PACK SIZE
2-CV	2mL, clear glass, crimp top, flat bottom, wide neck vial. With write-on patch.		500
2-RV	2mL snap cap vial with write on patch - clear.		500
09-FIV	900µL Micro+ crimp top, fused insert vial, clear.		500
06-PECV	600µL, polyethylene, crimp top vial with an internal tapered profile.		500
03-FIV	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
11-AC-ST101X	Aluminium crimp cap fitted with a pre-cut blue silicone/PTFE seal, with anti vacuum device.		500
11-PSN(B)-ST1X	11mm snap cap with pre-cut silicone/PTFE seal.		500

PART NUMBER	DESCRIPTION	6 x 32mm crimp top vial and cap	PACK SIZE
03-CVG	300µL, clear glass, crimp top, round bottom vial, made from inert Chromacol Gold grade glass, needs SV-S11A PTFE sleeve. See page 28.		500
8-AC-ST101X	Aluminium crimp cap fitted with a pre-cut blue silicone/PTFE seal.		500

PART NUMBER	DESCRIPTION	12 x 32mm wide neck screw top vials and cap	PACK SIZE
PEL	2mL screw top vial with a pre-cut white silicone/red PTFE seal.		100
2-SVW	2mL, clear glass, flat bottom, screw top, wide neck vial. With write-on patch.		500
09-FISV	900µL Micro+ screw top, fused insert vial, clear.		500
03-FISV	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
9-SC(B)-ST1X	9mm screw cap (blue) prefitted with pre-cut silicone/ PTFE seal.		500

PART NUMBER	DESCRIPTION	Wash vials and cap	PACK SIZE
10-CV	10mL, crimp top, headspace vial, used in Integral 4000.		125
6-CV	6mL, crimp top, headspace vial, used in the ISS100/200 and the Series 200.		125
20-PEPC5	20mm polyethylene snap cap.		250

