

All of these autosamplers were designed around Chromacol's own vial specifications. This guarantees a perfect fit every time when using Chromacol vials. The WebSeal system of 96 well titer plates and PTFE coated silicone elastomer closures are also compatible with Endurance and Reliance autosamplers. Instrument Select, pack SPL, includes 100 wide neck crimp top vials and the corresponding standard crimp caps.

Key to products ■ Vial ■ Cap ■ Seal ■ Combination Pack

Bold items in grey are our standard products, used in the combination packs.

PART NUMBER	DESCRIPTION	12 x 32mm crimp and snap cap vials and caps	PACK SIZE
SPL	2mL crimp top vial with a type 7 rubber/PTFE seal.		100
2-CV	2mL, clear glass, crimp top, flat bottom, wide neck vial. With write-on patch.		500
1.1-CTVG	1.1mL, clear glass, crimp top, tapered, wide neck vial. Made from inert Chromacol Gold grade glass. Needs TTS-312 PTFE support. See page 78.		500
09-FIV	900µL Micro+ crimp top fused insert vial, clear.		500
03-FIV	300µL Micro+, one piece, glass insert and clear vial combination. With write-on patch.		500
11-AC7	Aluminium crimp cap fitted with a type 7, rubber/PTFE seal. Shore hardness value 60. Other colours and seals available. See page 32.		500
2-RV	2mL snap cap vial with write on patch - clear.		500
11-PSN(B)-8RTI	11mm polyethylene snap cap pre-fitted with a red rubber/PTFE seal. Shore hardness value 58.		500

PART NUMBER	DESCRIPTION	8 x 40 and 6 x 32mm crimp top vial cap, seal	PACK SIZE
08-CV	800µL, clear glass, crimp top, flat bottom.		500
03-CVG	300µL, clear glass, crimp top, round bottom vial, made from inert Chromacol Gold grade glass, needs support sleeve SV-S3A. See page 28.		500
8-AC7	Aluminium crimp cap fitted with a type 7, rubber/PTFE seal. Shore hardness value 60. Other colours and seals available. See page 31.		1000

PART NUMBER	DESCRIPTION	12 x 32mm wide neck screw top vials and cap	PACK SIZE
2-SVW	2mL, clear glass, flat bottom, screw top, wide neck vial. With write-on patch.		500
09-FISV	900µL Micro+ screw top, fused insert vial, clear.		500
03-FISV	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
9-SC(B)-8RTI	Screw cap, blue, pre-fitted with type 8, rubber/PTFE seal. Shore hardness value 58. Other colours and seals available. See page 36.		500
1.1-STVG	1.1mL, clear glass, 8mm, screw top, tapered vial made from inert Chromacol Gold grade glass needs a TTS-312 support. See page 78.		500
8-SC-8RTI	8mm screw cap prefitted type 8 rubber/PTFE seal. Shore hardness value 57.		500

PART NUMBER	DESCRIPTION	22 x 38mm crimp top vials and cap	PACK SIZE
6-CV	6mL, clear glass, crimp top vial.		125
20-ACB	Aluminium crimp, no seal.		500
20-ST101	Blue silicone/PTFE, 1mm thick seal. Shore hardness value 30.		500

