

Chromacol's 12 x 32mm vials, both crimp and screw top versions are compatible with the above named autosamplers. A special PTFE sleeve, SV-TSP, has been produced to allow these units to operate with our 100, 200 and 300µL micro vials. **These sleeves however must not be used with the vial-heating unit.** Instrument Select product code TQL contains the recommended wide neck vials and caps.

Key to products ■ Vial ■ Cap ■ Combination Pack

Bold items in grey are our standard products, used in the combination packs.

PART NUMBER	DESCRIPTION	12 x 32mm crimp and snap cap vials and caps	PACK SIZE
TQL	2mL crimp top vial with a white silicone/red PTFE seal.		100
2-CV	2mL, clear glass, crimp top, flat bottom, wide neck vial. With write-on patch.		500
1.1-CTVG	1.1mL, clear glass, crimp top, tapered, wide neck vial. Made from inert Chromacol Gold grade glass. Needs TTS-312 PTFE support. <i>See page 78.</i>		500
09-FIV	900µL Micro+™ crimp top, fused insert vial, clear.		500
03-FIV	300µL, one piece, glass insert and clear vial combination. With write-on patch.		500
11-AC-ST15	Aluminium crimp cap fitted with a white silicone/red PTFE seal. Shore hardness value 57.		500
2-RV	2mL snap cap vial with write on patch - clear.		500
11-PSN(B)-STI	11mm polyethylene snap cap pre-fitted with a silicone/PTFE seal. Shore hardness value 57.		500

PART NUMBER	DESCRIPTION	12 x 32mm narrow and wide neck screw top vials and caps	PACK SIZE
2-SV	2mL, clear glass, flat bottom, screw top vial. With write-on patch.		500
1.1-STVG	1.1mL, clear glass, 8mm, screw top, tapered vial made from inert Chromacol Gold grade glass, needs a TTS-312 support. <i>See page 78.</i> Not suitable for older model Spectra-Physics autosamplers.		500
8-SC-ST15	Screw cap, with silicone/PTFE seal. Shore hardness value 57.		500
2-SVW	2mL, clear glass, flat bottom, screw top, wide neck vial. With write-on patch.		500
09-FISV	900µL Micro+ screw top, fused insert vial, clear.		500
06-PESV	600µL, polyethylene, 8mm, screw top vial with an internal tapered profile.		500
03-FISV	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
9-SC(B)-STI	Screw cap, blue, pre-fitted with silicone/PTFE seal. Shore hardness value 57. Other colours and seals available. <i>See page 36.</i>		500

PART NUMBER	DESCRIPTION	6 x 32mm crimp top vial and cap	PACK SIZE
03-CVG	300µL, clear glass, crimp top, round bottom vial, made from inert Chromacol Gold grade glass.		500
8-AC-ST15	Aluminium crimp cap fitted with a white silicone/red PTFE seal. Shore hardness value 57.		500

