

Varian GC Autosamplers 8000/8100/8410

All our 12 x 32mm standard and Micro+ vials are compatible with Varian autosamplers. We also have a wide range of vials suitable for microsampling and concentration by evaporation. Product code VAG contains the recommended screw top vials and the corresponding caps. For the full range of our GC Septa, please request a copy of our GC Septa brochure or see page 56.

Key to products ■ Vial ■ Cap ■ Combination Pack

Bold items in grey are our standard products, used in the combination packs.

PART NUMBER	DESCRIPTION	12 x 32mm crimp and snap cap vials and caps	PACK SIZE
2-CV	2mL, clear glass, crimp top, flat bottom, wide neck vial. With write-on patch.		500
1.1-CTVG	1.1mL, clear glass, crimp top, tapered vial made from inert Chromacol Gold grade glass, needs a TTS-314 PTFE support sleeve. See page 78.		500
09-FIV	900µL Micro+ crimp top, fused insert vial, clear.		500
06-PECV	600µL, polyethylene, crimp top vial with an internal tapered profile.		500
03-FIV	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
11-AC-ST15	Aluminium crimp cap fitted with a white silicone/red PTFE seal. Shore hardness value 57.		500
2-RV	2mL snap cap vial with write on patch - clear.		500
11-PSN(B)-STI	11mm polyethylene snap cap pre-fitted with a silicone/PTFE seal. Shore hardness value 57.		500

PART NUMBER	DESCRIPTION	12 x 32mm wide neck screw top vials and caps	PACK SIZE
VAG	2mL screw top vial with a white silicone/red PTFE seal.		100
2-SVW	2mL, clear glass, flat bottom, screw top, wide neck vial. With write-on patch.		500
09-FISV	900µL Micro+ screw top fused insert vial, clear.		500
03-FISV	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
9-SC(B)-ST1	9mm, blue, screw cap pre-fitted with a silicone/PTFE seal. Shore hardness value 57.		500
1.1-STVG	1.1mL, clear glass, screw top, tapered vial made from inert Chromacol Gold grade glass.		500
8-SC-ST15	8mm, screw cap pre-fitted silicone/PTFE seal. Shore hardness value 57.		500

PART NUMBER	DESCRIPTION	6 x 32mm crimp top vial and cap	PACK SIZE
03-CVG	300µL, clear glass, crimp top, round bottom vial, made from inert Chromacol Gold grade glass, needs SV-S3A support. See page 28.		500
8-AC-ST15	Aluminium crimp cap pre-fitted with a white silicone/red PTFE seal. Shore hardness value 57.		500

Chromseal GC Septa	Model	Septa size		PART NUMBER			Pack Size
		mm	ins.	High Temp	Long Life	Economy	
Varian	All packed column models	9.5	3/8	HT-9.5	LL-9.5	ECO-9.5	25
	3300,3400,3500,3600,	11	7/16	HT-11	LL-11	ECO-11	25
	3700,Vista	11	7/16	HT-11	LL-11	ECO-11	25

