

# Varian LC Autosamplers 9090/9095/PROSTAR

Our extended Micro+ range, with volumes of 300µL and 900µL are ideal for Varian autosamplers. Our ST-15 and ST-14 seals have the required degree of shore hardness for these instruments. There are many Chromacol vials that will operate with the Varian autosamplers, the most popular of which being our 12 x 32mm 2-CV & 2-SVW.

**Key to products** ■ Vial ■ Cap ■ Combination Pack

**Bold items in grey are our standard products, used in the combination packs.**

PART NUMBER	DESCRIPTION	12 x 32mm crimp and snap cap vials and caps	PACK SIZE
<b>2-CV</b>	2mL, clear glass, crimp top, flat bottom, wide neck vial. With write-on patch.		500
<b>1.1-CTVG</b>	1.1mL, clear glass, crimp top, tapered vial made from inert Chromacol Gold grade glass, needs a TTS-314 PTFE support sleeve. <a href="#">See page 78.</a>		500
<b>09-FIV</b>	900µL Micro+ crimp top, fused insert vial, clear.		500
<b>06-PECV</b>	600µL, polyethylene, crimp top vial with an internal tapered profile.		500
<b>03-FIV</b>	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
<b>11-AC7</b>	Aluminium crimp cap fitted with a type 7, rubber/PTFE seal. Shore hardness value 60.		500
<b>2-RV</b>	2mL snap cap vial with write on patch - clear.		500
<b>11-PSN(B)-8RT1</b>	11mm polyethylene snap cap pre-fitted with a red rubber/PTFE seal. Shore hardness value 58. Other colours and seals available. <a href="#">See page 25.</a>		500

PART NUMBER	DESCRIPTION	12 x 32mm wide neck screw top vials and caps	PACK SIZE
<b>VAL</b>	2mL screw top vial with a white silicone/red PTFE seal.		100
<b>2-SVW</b>	<b>2mL, clear glass, flat bottom, screw top, wide neck vial. With write-on patch.</b>		<b>500</b>
<b>09-FISV</b>	900µL Micro+ screw top, fused insert vial, clear.		500
<b>03-FISV</b>	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
<b>9-SC(B)-ST1</b>	<b>Screw cap, blue, pre-fitted with silicone/PTFE seal. Shore hardness value 57.</b> <b>Other colours and seals available. <a href="#">See page 36.</a></b>		<b>500</b>
<b>1.1-STVG</b>	1.1mL, clear glass, 8mm, screw top, tapered vial made from inert Chromacol Gold grade glass. Needs a TTS-312 support. <a href="#">See page 78.</a>		500
<b>8-SC-ST15</b>	8mm, screw cap pre-fitted silicone/PTFE seal. Shore hardness value 57.		500

PART NUMBER	DESCRIPTION	6 x 32mm crimp top vial and cap	PACK SIZE
<b>03-CVG</b>	300µL, clear glass, crimp top, round bottom vial, made from inert Chromacol Gold grade glass, needs an SV-S3A support. <a href="#">See page 28.</a>		500
<b>8-AC7</b>	Aluminium crimp cap fitted with a type 7, rubber/PTFE seal. Shore hardness value 60. Other colours and seals available. <a href="#">See page 29.</a>		1000

