

Waters 48/717/717 Plus

The Chromacol 4-SV(A), 15 x 45mm, 4mL screw top, amber vial is also an ideal vial for this autosampler. When using our extensive range of PTFE sleeves, many of our micro vials may be used with these autosamplers.

Our Sci-Vi system vials, including the 01-CVG vial, may be used as reduced volume inserts in a 4-SV vial with an S-50 spring. 07-CPV(A) vial can be used as a reduced volume amber insert with an S-15 spring. Our preferred solution is the SV-S15, a new one piece low cost plastic support which replaces the above mentioned 4mL vial. This easy and simple to use support uses any one of the Sci-Vi system vials, i.e. 03-CVG, 02-CTVG, 02-CTV(A) and the 01-CVG.

Key to products ■ Vial ■ Cap

Bold items in grey are our standard products

PART NUMBER	DESCRIPTION	15 x 46mm screw top vial and cap	PACK SIZE
4-SV	4mL, clear glass, screw top, flat bottom vial.		500
12-SC-ST2	Screw cap, black, pre-fitted with a 2mm thick, silicone/PTFE seal.		500
3.5-HRSV	3.5mL screw top vial - High Recovery - uses 13mm screw caps.		500
13-SC-ST15	13mm screw cap prefitted silicone/PTFE seal.		500

PART NUMBER	DESCRIPTION	12 x 32mm crimp top vials and cap	PACK SIZE
4-CV	4mL, clear glass, crimp top, flat bottom vial, 15 x 46mm.		500
2-CV	2mL, clear glass, wide neck, crimp top, flat bottom vial. With write-on patch.		500
09-FIV	900µL Micro+ crimp top, fixed insert vial, clear.		500
03-FIV	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
11-AC-ST15	Aluminium crimp cap fitted with a white silicone/red PTFE seal. Shore hardness value 57.		500

PART NUMBER	DESCRIPTION	6 x 32mm crimp top vial and cap	PACK SIZE
03-CVG	300µL, clear glass, crimp top, round bottom vial, made from inert Chromacol Gold grade glass, needs support sleeves SV-S11A and WS-2 or SV-S15. See page 28.		500
8-AC-ST15	Aluminium crimp cap fitted with a white silicone/red PTFE seal. Shore hardness value 57.		500

Note

Standard 12 x 32mm vials require the use of the WS-7 support sleeve (see page 78), when used on the 717 and 717 Plus.

