

# Waters 96/Alliance/Acquity

The 96 position Waters autosampler does not use standard 12 x 32 mm vials. The Chromacol I.2-CWV and I-CWV vials are specially designed to meet the requirements of this instrument. When using small samples we would recommend the use of the I-CWV vial because of its tapered base.

**Key to products** ■ Vial ■ Cap ■ Combination Pack

**Bold items in grey are our standard products, used in the combination packs.**

## Waters 96

PART NUMBER	DESCRIPTION	8 x 40mm crimp top vials and cap	PACK SIZE
<b>I.2-CWV</b>	<b>1.2mL, clear glass, crimp top, flat bottom vial.</b>		<b>500</b>
<b>I-CWV</b>	1mL, clear glass, crimp top, tapered vial.		500
<b>8-AC-ST15</b>	<b>Aluminium crimp cap fitted with a white silicone/red PTFE seal. Shore hardness value 57.</b>		<b>500</b>

PART NUMBER	DESCRIPTION	8 x 40mm neckless vials and caps	PACK SIZE
<b>I-NWV</b>	1mL, clear glass shell vial. Uses 8-NPWP polyethylene plug closure.		500
<b>I-NWV-C</b>	1mL neckless clear glass vial, with polyethylene caps.		200
<b>I-NWV(A)-C</b>	1mL neckless amber glass vial, with polyethylene caps.		200
<b>I-PPNWV-C</b>	1mL PP neckless vial for Waters 96 with a polyethylene cap.		250

## Alliance/Acquity

For these autosamplers the Instrument Select code WAL includes our wide neck 2-SVW vials and our 9-SC(B)-STI silicone/PTFE seals. Waters now recommend the use of a bonded cap with some Alliance/Acquity instruments.

PART NUMBER	DESCRIPTION	12 x 32mm crimp and snap cap vials and caps	PACK SIZE
<b>2-CV</b>	2mL, clear glass, crimp top, flat bottom, wide neck vial. With write-on patch.		500
<b>1.1-CTVG</b>	1.1mL, clear glass, crimp top, tapered vial made from inert Chromacol Gold grade glass, needs a TTS-314 PTFE support sleeve. <a href="#">See page 78.</a>		500
<b>09-FIV</b>	900µL Micro+ crimp top, fused insert vial, clear.		500
<b>06-PECV</b>	600µL, polyethylene, crimp top vial with an internal tapered profile.		500
<b>03-FIV</b>	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
<b>11-AC-ST15</b>	Aluminium crimp cap fitted with a white silicone/red PTFE seal. Shore hardness value 57.		500
<b>2-RV</b>	2mL snap cap vial with write on patch - clear.		500
<b>11-PSN(B)-STI</b>	11mm polyethylene snap cap pre-fitted with a silicone/PTFE seal. Shore hardness value 57.		500

PART NUMBER	DESCRIPTION	12 x 32mm wide neck screw top vials and cap	PACK SIZE
<b>WALB</b>	2mL screw top vial with a bonded white silicone/red PTFE seal.		100
<b>2-SVW</b>	<b>2mL, clear glass, flat bottom, screw top, wide neck vial. With write-on patch.</b>		<b>500</b>
<b>09-FISV</b>	900µL Micro+ screw top, fused insert vial, clear.		500
<b>03-FISV</b>	300µL Micro+, one piece, clear glass insert and vial combination. With write-on patch.		500
<b>9-SC(B)-STI</b>	Screw cap, blue, pre-fitted with silicone/PTFE seal. Shore hardness value 57. Other colours and seals available. <a href="#">See page 36.</a>		500
<b>9-SC(BLK)-BSTI</b>	<b>9mm screw cap (black) bonded silicone/PTFE liner. Recommended by Waters.</b>		<b>500</b>
<b>9-SC(GY)-BSTIX</b>	9mm screw cap (grey) bonded silicone/PTFE liner pre-cut. Recommended by Waters.		500