

## Chromacol Chromatography Consumables

Building on our history of delivering innovative solutions, we have a number of new products to help you improve your chromatography experience.

### 11mm Magnetic Crimp Caps



These soft crimp magnetic caps allow the use of 2mL crimp vials for sample pre-treatment and preparation on autosamplers requiring magnetic transport such as the CTC CombiPAL.

See page 32



### New High Recovery Storage Vials

Chromacol's new range of high recovery storage vials contain a tapered reservoir. This provides larger sample capacities whilst giving maximum sample recovery and prevents any waste of precious samples.

See page 42

### Storage Vials

For general storage applications these dram vials take caps with standard GPI threads. Both injection caps with piercable seals or solid PTFE lined storage caps may be supplied separately, and also as part of convenience packs.



See page 69

### Bonded Caps

Chromacol bonded caps are designed to fit the full range of Chromacol 9mm screw top vials. Using a reagent free bonding process, the seal is bonded into the polypropylene cap. They are available in either a standard or pre-cut version.



See page 36

## Thermo Scientific Chromatography Columns and Consumables

We now offer a number of Thermo Scientific accessories for chromatography and mass spectrometry.

### Thermo Scientific HyperSep SPE Columns



HyperSep™ products have been developed for rapid effective and economical sample preparation.

See page 52

### Thermo Scientific HPLC Columns



Chromacol now supplies columns packed with Thermo Scientific Hypersil™ materials. This includes the Hypersil GOLD™ range of base deactivated phases, including 1.9µm diameter particle for ultra high pressure applications.

See page 44